

REEP2-receptor accessory protein 2 Gene

Reference: PTXBC006218

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C5orf19; SGC32445; SPG72; Yip2d; receptor expression-enhancing protein 2; receptor accessory protein 2		

For research use only.

GENE DESCRIPTION	receptor accessory protein 2
GENE ID	51308
ORIGIN SPECIES	Human
ACCESSIONS	BC006218
VECTOR	pDONR223
NCBI SYMBOL	REEP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgtctctggtatctctcgcctggtgtgtcatcttggcacctgtaccagcctattctctacaaggccgtgaagacaaaaacgtgaaggaaatgtgaaatggatgatgtactggatcgtcttgccttctcaccacggccgagacgtcacgatatagtgctctcgtgttcc
ccttctacttgaactgaagatgccttctgtgataggtgtgtgtcccttacccaagggtccagcgtgtgtctaccgcaagttcgtgcaccaacgtgttccacaaggagaaggagatcagagtagcatcacgcaggcccagacaaagctatgagaccatgatgaggtgggca
agaggggctgaaccttgcgcgaatgctgcagtcacagctgccccaagggtgtgtgtcagagaagctccgcagcttcagcatgcaggacctgacctgatccgggacgaggacgacagctccccctgagaggcctgacggccgctccgacccagccctggcagcctctgga
caccatcaggacttaggagatgacctgacctgagcttaaggtccagcaaacccggcagattcccgacagaggcttctgaggatgacatgggagacaaagctcccaaggggccaaacccatcaaaaagcgcccaagctgagccactggcttccaagacactgaagacc
ggcccaagaagaagaccttgcgggggagcagcttga