

ZNF839-zinc finger protein 839 Gene

Reference: PTXBC006222

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf131; zinc finger protein 839; renal carcinoma antigen NY-REN-50		
For research use only.			

GENE DESCRIPTION	zinc finger protein 839
GENE ID	55778
ORIGIN SPECIES	Human
ACCESSIONS	BC006222
VECTOR	pDONR223
NCBI SYMBOL	ZNF839
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggttaagaaatgtccaagattacctcagtagttctgtgtgtccaggagaccctggaaataacaatgataagggtgctgagtcattaggaatcacagaattctacggaagaagaataacacccagacaacctggaccaagcacctcagccgagacatggatggggagcagctagagggagctagcagcgagaagagggaaactgagggctgaggagggagctggcctcagtgaaaaggccagagaagagccctgtccaacgataccactgaatctctgtgccaacagcagagggcgggagaagccagggcccttgcatttggccgtgtgacaatagtgtctcaggaggaggacattgtcacagtgcagtgatgcagaggggctgctgcgatgggcccgtag