

NUDT16L1-nudix (nucleoside diphosphate linked moiety X)-type motif 16-like 1 Gene

Reference: PTXBC006223

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SDOS; protein syndesmos; NUDT16-like protein 1; nudix (nucleoside diphosphate linked moiety X)-type motif 16-like 1; syndesmos; nudix hydrolase 16 like 1		

For research use only.

GENE DESCRIPTION	nudix (nucleoside diphosphate linked moiety X)-type motif 16-like 1
GENE ID	84309
ORIGIN SPECIES	Human
ACCESSIONS	BC006223
VECTOR	pDONR223
NCBI SYMBOL	NUDT16L1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgacggcggttcggagctgaagcagatcagccgggtggaggcgaatgcgcctagggcggtggagccactcgtgccagccatgctgtacgcccaaccctggcgagctcttcggccgcatccccatgcgttctcggtgctgatcagatcggttcgacgggctg
ctgggcttccccggggcttcgtggaccggcgttctggtcgtggaggacggcctgaaccgggtgctgggctgggctgggctggcctgcgcctcaccgagggcgactacctgagctcgacctgaccgagggccacaccgctcgtggcgacctgtacgcggcagctgac
gctggagcagctgcacgccgtggagatcagcgggtgactcgcgcaccacggcctggaggtgctgggctcgtggggtcccgctgtacaccagaaggaccgagtcggaggctccccaaacttctgagcaacgccttcgtgagcacggctaagtccagctcctcttggccctc
aaggtgctcaacatgatcccgaggagaagctggttgaggccctgctgcagccaccgagaagcagaagaaggccctggagaagttgctcccgccctcttga