



Reference: PTXBC006232

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDC45L; CDC45L2; MGORS7; PORC-PI-1; cell division control protein 45 homolog; CDC45 cell division cycle 45 homolog; CDC45 cell division cycle 45-like; CDC45-related protein; cell division cycle 45 homolog; cell division cycle 45-like 2; human CDC45; cell		

GENE DESCRIPTION	CDC45 cell division cycle 45-like (<i>S. cerevisiae</i>)
GENE ID	8318
ORIGIN SPECIES	Human
ACCESSIONS	BC006232
VECTOR	pDONR223
NCBI SYMBOL	CDC45L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]