

## RHBDD2-rhomboid domain containing 2 Gene

Reference: PTXBC006234

|                |                                                                                                              |                     |                |
|----------------|--------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                          | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                         | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | NPD007; RHBDL7; rhomboid domain-containing protein 2; rhomboid, veinlet-like 7; rhomboid domain containing 2 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | rhomboid domain containing 2                                                                      |
| <b>GENE ID</b>          | 57414                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC006234                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | RHBDD2                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgctgggagtcaccaccgtccgttctcggatgagcgggccctggtgtttgcatggtgtgccccagtcctggttccgtggctcctgctgggtgctcgtgggtcattccccagaccttttctcagtaattgtctgctgggtgtccatcgggtggcctatggcctcacctactgtattcc  
atcgacctctcagagcgagtggtcgtgaagctcgaatcagaccttccccctcagcctgatgaggagatgccgtgttcaagtacgtctcaggggtctcagccgagagggagggcagccagagccggaaactgaaccgggtgctgtctacccccacagagctgccacctcacct  
gtcccaagccaccctgtgtccagacgcagcagccagtggtcagaagctggcctcctggccctcctgcacccccgggcacatgccaccttgctcctgaccagcctgctccgctgtgtatgtgcagaaccatttgggtccaaacccacctctccagtgcttaccagcttct  
gcgggcacctccctggcctcagccccacgctgtgaacagccctggcacggtgtattctggggccttgggcacaccaggggctgcaggctccaaggagtcctcagggtccccatgccctga