



Reference: PTXBC006242

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	major facilitator superfamily domain-containing protein 9; major facilitator superfamily domain containing 9		

For research use only.

GENE DESCRIPTION	major facilitator superfamily domain containing 9
GENE ID	84804
ORIGIN SPECIES	Human
ACCESSIONS	BC006242
VECTOR	pDONR223
NCBI SYMBOL	MFSD9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1