

RGS16-regulator of G-protein signaling 16 Gene

Reference: PTXBC006243

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	A28-RGS14; A28-RGS14P; RGS-R; regulator of G-protein signaling 16; hRGS-r; regulator of G-protein signalling 16; retinal-specific RGS; retinally abundant regulator of G-protein signaling		

For research use only.

GENE DESCRIPTION	regulator of G-protein signaling 16
GENE ID	6004
ORIGIN SPECIES	Human
ACCESSIONS	BC006243
VECTOR	pDONR223
NCBI SYMBOL	RGS16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgccgcaccctggccgcttccccaccacctgcctggagagagccaaagattcaagacacgtctggggatcttttcacaaatcagagctgggctgcgatactgggagtactggcaagttcagtggtggcgagtaaacacagcaaaagagaatagaacttctcagaagatgtctg
gggtggagagagctgttcgacctgctgtagcagtaaaatggagtggtgcttccacgcttctgaagacagagttcagtgaggagaacctggagttctgctggcctgtgaggagttcaagaagatccgatcagctaccaagctggcctccagggcacaccagatctttgaggagt
tcatttcagtgaggccctaaagaggtaacattgacctgagaccccgagctgacgaggaatgaacctgcagactgccacagccacatgcttgcggtcaggggaagacacgtacacctgatggagaaggactcctaccacgcttctgaagtcgctgcttaccgggacctg
gctgcccaagcctcagccgctctgccactctgtccagctgcagcctggacgagccctcacacacctga