

C15orf41-chromosome 15 open reading frame 41 Gene

Reference: PTXBC006254

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C15orf41; HH114; chromosome 15 open reading frame 41		
For research use only.			

GENE DESCRIPTION	chromosome 15 open reading frame 41
GENE ID	84529
ORIGIN SPECIES	Human
ACCESSIONS	BC006254
VECTOR	pDONR223
NCBI SYMBOL	C15orf41
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcggcttatactggagaggtttctacaggaacacagaggaactccacctccaagtctattataatagtagtctacgggaccctctcagattccagatggagttctagcaaatcaggctctatcagtcattgtgaacgactgctgttacggaccactagtggaactgcatcaagcatgc
cattggctcatgagcatgaggtcctgctgagagactgctctagagaaaaacctgtcctcctagatgaagatcagcttcgtgcaaaagggtatgacaaaacaccagacttcatttacaagtaccagttgctgtagaagggcacataattcactggattgaaagcaaacctcatttggtgatga
atgtagccaccacgcctacctgcatgaccagttctggagctactggaatagatttggccaggcttagtcatctattggtatggattatccaggagctggactgcaaccgggaaaggggcatcctgctcaaaagcctgtttcccccagaacattgtcaccttatgccacagcatagcttga