

VPS25-vacuolar protein sorting 25 homolog (S. cerevisiae) Gene

Reference: PTXBC006282

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ESCRT-II complex subunit VPS25; DERP9; FAP20; vacuolar protein-sorting-associated protein 25; ELL-associated protein of 20 kDa; dermal papilla-derived protein 9; vacuolar protein sorting 25 homolog		

For research use only.

GENE DESCRIPTION	vacuolar protein sorting 25 homolog (S. cerevisiae)
GENE ID	84313
ORIGIN SPECIES	Human
ACCESSIONS	BC006282
VECTOR	pDONR223
NCBI SYMBOL	VPS25
LENGTH	Proteogenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgatgagttcgagtgccgtggcagatcgtccacccttttacgttaacaaccgaatgtggacactggcagaagcagctggccgctggtgctcgtgtctctctcttgcgcctgcacaacagtcacgatgacggatggaagctcaggagagcccgtcttcaa
caacgtcaagctacagcgaagcttctgtgagtcgatccagattgtattagaggaaactgaggaaagaaaggaaactcagtggttggataagagcaagtcagcttctgatcatgtggcggaggccagaagaatggggaaactatctatcagtttccaggagtgccagaa
caactcgtctttaccctgtatgaactgactaatgggaagacacagaggatgaggattccacggcctggatgaagccactctactcgggctctcaggccctacagcaggagcacaaggccgagatcatcactgtcagcagtgccgagcgtcaagttctctag