



Reference: PTXBC006283

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SDR47C1; retSDR3; 17-beta-hydroxysteroid dehydrogenase 14; 17-beta-HSD 14; 17-beta-hydroxysteroid dehydrogenase DHRS10; dehydrogenase/reductase (SDR family) member 10; retinal short-chain dehydrogenase/reductase 3; retinal short-chain dehydrogenase/reduct		

For research use only.

GENE DESCRIPTION	hydroxysteroid (17-beta) dehydrogenase 14
GENE ID	51171
ORIGIN SPECIES	Human
ACCESSIONS	BC006283
VECTOR	pDONR223
NCBI SYMBOL	HSD17B14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1