

AXIN2-axin 2 Gene

Reference: PTXBC006295

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AXIL; ODCRCS; axin-2; axin-like protein; axis inhibition protein 2; conductin; axin 2		
For research use only.			

GENE DESCRIPTION	axin 2
GENE ID	8313
ORIGIN SPECIES	Human
ACCESSIONS	BC006295
VECTOR	pDONR223
NCBI SYMBOL	AXIN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtagcgtatgttggtgacttgcctccggacccagcagcagcttcgtgaggatgccccgccccccagtgccaggggaagaaggggagacccaccgtgtcagccaggggtgggcaagggccaggtcaccaaacccatgcctgtcttccaacaccaggcggaacga
agatgggttggggagccggagggggcgccatcctcggtatccctctgacccgggtggaccaagtccttacactccttattggcgcatcaagacggtgcttacctgtccgaacttctctggagaggagaaaatgcgtggatagcttagacttctgtttgcctgcaatggatcaggcagat
gaacctgaaggataccaaaactttacgagtagccaatgactttgtcatgacattttgttctacttatactgtaaaattatgcattataaagagttcatttaaggaaaattacttggtaataaa