

C1orf170-chromosome 1 open reading frame 170 Gene

Reference: PTXBC006300

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf170; PGC-1 and ERR-induced regulator in muscle protein 1; PGC-1 and ERR-induced regulator in muscle 1; PPARGC1 and ESRR-induced regulator in muscle 1; peroxisome proliferator-activated receptor gamma coactivator 1 and estrogen-related receptor-induce		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 170
GENE ID	84808
ORIGIN SPECIES	Human
ACCESSIONS	BC006300
VECTOR	pDONR223
NCBI SYMBOL	C1orf170
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgcctggtgttagcttttgcacactgggctgtgaggacgtcagatccgcatacccccagacgctggaaaacagccttgcctggccaacgtcggcaccatctctgccatccgtacttccgccggcaggcgggggcaaggggccgcagccacagccccagccctcctag