

MSC-musculin (activated B-cell factor-1) Gene

Reference: PTXBC006313

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ABF-1; ABF1; MYOR; bHLHa22; musculin; activated B-cell factor 1, homolog of mouse musculin; activated B-cell factor-1; class A basic helix-loop-helix protein 22		

For research use only.

GENE DESCRIPTION	musculin (activated B-cell factor-1)
GENE ID	9242
ORIGIN SPECIES	Human
ACCESSIONS	BC006313
VECTOR	pDONR223
NCBI SYMBOL	MSC
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccacgggctcgtgagtgatccggaggagatggagctcgggggctgcagcgggagctacccgggtcccgctcccaagagggccgcccctccgctgtagagcgcagctacgcctcggcagtgacaactcgtcggcagaggaggagacccgacggcgaggaggagc
gctgcgtctgggcacagccggcagcgcgggaaggctgcaagaggaagcggcccgctgtggtgctggggcgggcgagctgtagcgcggcggtggtggcaagaagccctccggccaagggtcagccgcagagtgcaagcagtcgcagcggaacgcggccaacgcc
cgtgagcgtgccggatgctgctgagcaagccttctccaggctcaagaccagcctgcccgtgggtgcccccgacactaagctctccaagctggacagctccggctggcttcagttacatcgtcactcgcggcagctgtgtaggagaccgtatgagaacggctacgtgca
ccagtgaaacctgacatggccattcgtgctcgggaagaccggactctgacacaaagaagtttccgcagccaacagactatgtggaaccacgcttaa