

AP2S1-adaptor-related protein complex 2, sigma 1 subunit Gene

Reference: PTXBC006337

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AP17; CLAPS2; FBH3; FBHOk; HHC3; AP-2 complex subunit sigma; HA2 17 kDa subunit; adaptor protein complex AP-2 subunit sigma; clathrin assembly protein 2 sigma small chain; clathrin coat assembly protein AP17; clathrin coat-associated protein AP17; clathri		

For research use only.

GENE DESCRIPTION	adaptor-related protein complex 2, sigma 1 subunit
GENE ID	1175
ORIGIN SPECIES	Human
ACCESSIONS	BC006337
VECTOR	pDONR223
NCBI SYMBOL	AP2S1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatccgctttatctcaccagaacggcaggaagacgcgcctggccaagtgtacatgcagtttgatgatgagaaacagaagctgacgaggaggtgcatgccgtggtcaccgtccgagacgccaacacaccaactttgtggagttccggaactttaagatcatttacgcc
gctatgctggcctctactctgcatctgtgtggaatgcaatgacaacacctggcttacctggaggccattcacaactctgtggaggtcttaaacgaatattccacaatgtctgtgaactggacctggtgttcaacttctacaaggtttacacggctgtggagagatgttctgctggcgaatc
cgagagaccagccagacgaaggtgctgaacagctgctgatgctacagtccttgagtgga