

## NCK1-NCK adaptor protein 1 Gene

Reference: PTXBC006403

|                |                                                                                                                                                                                         |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                     | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                    | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | cytoplasmic protein NCK1; NCKalpha; nck-1; NCK tyrosine kinase; SH2/SH3 adaptor protein NCK-alpha; melanoma NCK protein; non-catalytic region of tyrosine kinase; NCK adaptor protein 1 |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | NCK adaptor protein 1                                                                             |
| <b>GENE ID</b>          | 4690                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC006403                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | NCK1                                                                                              |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggcagaagaagtggtgtagtgcacaaattgattatgtgcccacaagaacaagagtggtgacatcaagaagaatgagagattatgcttctggtgagtggttcgaaattccatgaataaacaggtttgtgcttctaactatgtggaaggaaaaacagtgc  
tcggaaagcatctatttgaaaaacctaaggatccttaggcattggaaaagtgaagaagaacctaagtggtgacatctctcctgctgatagattttgtgacccagggaacgtctctatgacctcaacatgccgcttatgtgaatttaactacatggctgagagagagatga  
attatcattgataaaggggacaaaggtgatcgtcatggagaatgcagtgatgggtggtggtggtgtagctacaatggacaagtggtggttccctcaactatgtaactgaagaaggtgacagtcctttgggtgacctgtgggtctctgtagagaaattagcagcagtcgtaatac  
ctaaatactggcgaagtgttgcatgtgttacaggctctttacccattcagctcatctaatgatgaagaacttaattcgagaaggagatgtaattgattgaaacacctgaaatgacccagagtggtggaatgcaggaaatcaatggatgtgttggtctgtaccataaactatgta  
ccgttatgcagaataatccattaactcaggtttggaaccatcacctccacagtggtgattacattaggccctcactcactggaagtttgcggcaatccttggtattatggcaagtcaccagcagcaagcagaatggcattaaatgaaggagacatgaagggttctcctcattcgtgatag  
tgaatcttcgcaaatgatttctcagtcactaaagcacaagggaagaacacagcattttaaagtcacactaaagagagactgtctactgcatggcgagcgttaattcagcaccatggagaacttgtagaacattacaaaaggcaccatttttacaagtgaacaaggagaaaaattatct  
tgtcaagcatttatcatga