



Reference: PTXBC006411

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CANDF9; IL17-RL; IL17RL; interleukin-17 receptor C; IL-17 receptor C; IL-17RL; IL17F receptor; IL17Rhom; interleukin-17 receptor homolog; interleukin-17 receptor-like protein; zcytoR14; interleukin 17 receptor C		

For research use only.

GENE DESCRIPTION	interleukin 17 receptor C
GENE ID	84818
ORIGIN SPECIES	Human
ACCESSIONS	BC006411
VECTOR	pDONR223
NCBI SYMBOL	IL17RC
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr