

CCDC115-coiled-coil domain containing 115 Gene

Reference: PTXBC006429

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDG2O; ccp1; coiled-coil domain-containing protein 115; coiled-coil domain containing 115		

For research use only.

GENE DESCRIPTION	coiled-coil domain containing 115
GENE ID	84317
ORIGIN SPECIES	Human
ACCESSIONS	BC006429
VECTOR	pDONR223
NCBI SYMBOL	CCDC115
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccaagtggtagggccctgcagatgcttcccatggagccccaggtctgcctccacgccagcaggccccagggagctcagaagttcaagtggtgagagctggtgtccacgccagaggaggtggggcctcgcgaagcaggtctgcggaggcgcaaggc
cccactaagacccagaaaccggagtcctctgagggccctcaggacccctgaactggtttggaatcctagtctcacaagtcacgtcaggtcgaagcaagcttcgggatggcctgcagctggccgcagacatagccagcctccagaaccgattgactgggtcgaagccagctccg
gggactccaagagaaactcaagcagctggagcctggggctgctga