

PCDHGC3-protocadherin gamma subfamily C, 3 Gene

Reference: PTXBC006439

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PC43; PCDH-GAMMA-C3; PCDH2; protocadherin gamma-C3; cadherin-like 2; protocadherin 43; protocadherin-2; protocadherin gamma subfamily C, 3		

For research use only.

GENE DESCRIPTION	protocadherin gamma subfamily C, 3
GENE ID	5098
ORIGIN SPECIES	Human
ACCESSIONS	BC006439
VECTOR	pDONR223
NCBI SYMBOL	PCDHGC3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtcccagagcctggaggagcggactgcaagccccgcccaacacggactggcggttctctcagggccagagaccggcaccagcggtcccaaatggcgatgacacggcacctggcccaacaaccagtttgacacagagatgctgcaagccatgatcttggcgccgag
tgaagctgctgaggagctccacctgggaggggtgcccgcacatgggattgagcgccgctacggacccagttcacctgcagcacgtgcccgactaccgccagaatgtctacatccaggcagcaatgccacactgaccaacgcagctggcaagcgggatggcaaggcc
ccagcaggtggcaatggcaacaagaagaagtcgggcaagaaggagaagaagtaa