

CKS2-CDC28 protein kinase regulatory subunit 2 Gene

Reference: PTXBC006458

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CKSHS2; cyclin-dependent kinases regulatory subunit 2; CDC28 protein kinase 2; CKS-2; CKS1(S. cerevisiae Cdc28/Cdc2 kinase subunit) homolog-2; CDC28 protein kinase regulatory subunit 2		

For research use only.

GENE DESCRIPTION	CDC28 protein kinase regulatory subunit 2
GENE ID	1164
ORIGIN SPECIES	Human
ACCESSIONS	BC006458
VECTOR	pDONR223
NCBI SYMBOL	CKS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccacaagcagatctactctcgacaagtactcgacgaacactacgagtagccgcatgttatgttaccagagaacttccaacaagtacctaatactatctgatgtctgaagaggagtgaggagacttggtgtccaacagagtctaggctgggttcattacatgattcatgagc
cagaaccacatattctctcttttagacgaccttccaaaagatcaacaaaatga