

## CSNK1E-casein kinase 1, epsilon Gene

Reference: PTXBC006490

|                |                                                                                                |                     |                |
|----------------|------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                            | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                           | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | CKIepsilon; HCKIE; casein kinase I isoform epsilon; CKI-epsilon; CKIε; casein kinase 1 epsilon |                     |                |

**For research use only.**

---

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | casein kinase 1, epsilon                                                                         |
| <b>GENE ID</b>          | 1454                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC006490                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | CSNK1E                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggagctacgtgtgggaacaagtaccgcctgggacggaagatcgggagcgggtcctcggagatatctacctgggtgccaatcgccctctggtgaggaagtcgccatcaagctggagtgtgtgaagacaaagcaccgccagctgcacatcgagagcaagttctacaagatgatc  
aggggtgctggggatcccgatcaagtggtgctggagctgagggcgactacaacgtgatgtgatgagctgtgtgggctgagcctcgaggacctgttcaactctgttccgcaaatcagcctcaagacgggtgctgttggccgaccagatgatcagccgatcagtagtatatcca  
ctccaagaacttcatccacgggacgtcaagcccgacaacttctcatgggctggggaagaaggcaacctgtgtctacatcatcgacttggcctggccaagaagtaacgggacgcccgcaccaccagcacattccctaccgggaaacaagaacctgacgggcacggccgcta  
cgcttccatcaacacgcacctgggcatggagcaaacgctgagatgacctggagagcctgggctacgtgtctatgtacttcaactgggctcctcgccctggcaggggctcaaaagcagcaccacagcgcagagaagatgatgaacggatcagcgagaagaagatgtcaacgccatcga  
ggctccttcaaaaggctatccctccgaattctcaacatacctcaacttctgctcctcctgggtttgacgacaagcccgactactcttactacgtcagctcttccgaacctcttccaccggcagggtcttctctatgactacgtcttggactggaacatgtgaaatcgggtgcagccggaat  
cccaggatgtgacgggagcggcgagaacacgaacgcgagagaggtggggcagctacgggggtccgcgacccgagccctgccccctggcccaccacggggccactgccaacgggtccgcagtgccggcagccgtggcttccacgccagcctcccgatccag  
ccggctggcaatacttctccagagcgatctcgcggtcgaccgggagagggaaggtgagtagggctgcacaggggtgcggccgcaacgtctcctcagacctcactggcgcggaagaggtctcccgatcccgctcagacagacaagtggtccatttgacatctcgggaagt  
ga