



Reference: PTXBC006495

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CML33; FARSL; FARSLA; FRSA; PheHA; phenylalanine--tRNA ligase alpha subunit; pheRS; phenylalanine tRNA ligase 1, alpha, cytoplasmic; phenylalanine--tRNA ligase alpha chain; phenylalanine-tRNA synthetase alpha-subunit; phenylalanine-tRNA synthetase-like, a		

For research use only.

GENE DESCRIPTION	phenylalanyl-tRNA synthetase, alpha subunit
GENE ID	2193
ORIGIN SPECIES	Human
ACCESSIONS	BC006495
VECTOR	pDONR223
NCBI SYMBOL	FARSA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr