

HOXA9-homeobox A9 Gene

Reference: PTXBC006537

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	homeodomain protein HOXA9; ABD-B; HOX1; HOX1.7; HOX1G; homeobox protein Hox-A9; homeobox protein Hox-1G; homeobox A9		

For research use only.

GENE DESCRIPTION	homeobox A9
GENE ID	3205
ORIGIN SPECIES	Human
ACCESSIONS	BC006537
VECTOR	pDONR223
NCBI SYMBOL	HOXA9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaccactggggccctgggcaactactacgtggactgttctgtggggcgccgacgcccggatgagctgagcgttggccgctatgcgcccgggaccctggccagcctccccggcaggcgccgacgctggccgagcaccgccagcttcagccgtgcagcttcagtcctccaa
ggcgacgggtttggcgccctgtggaacccagtgacgcggcgccgccaacgtgtaccgctgctgtgtaccaccaccatcaccaccaccctacgtgcacccccaggcgcccggtggcgccggcgccggacggcagggtacatgcctctggttgagccccacggccgg
tgcgctctcttcggggcttgcctccagccggccttatggcattaaactgaaccgctgtcggccagaaggggtgactgtccacgcttgacactcacactttgtccctgactgactatgcttgggttctctccagttgalagagaaaaaacaccagcgaaggcgcttctgaaaaac
aatgctgagaatgagagcgccgagacaagcccccatcgtcccaataaccagcagcgaactggcttcatgcgcgtccactcggaagaaagcggtgccctatacaaaacaccagaccctggaactggagaagagtttctgttcaacatgtacctaccaggaccgcaggtacg
agggtgctcgactgtcaacctcaccgagaggcagggtcaagatctgttcagaaccgcaggatgaaatgaagaaatcaacaagaccgagcgaagacgagtgga