

PTRH2-peptidyl-tRNA hydrolase 2 Gene

Reference: PTXBC006807

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BIT1; CFAP37; CGI-147; IMNEPD; PTH; PTH 2; PTH2; peptidyl-tRNA hydrolase 2, mitochondrial; bcl-2 inhibitor of transcription 1; cilia and flagella associated protein 37; peptidyl-tRNA hydrolase 2		

For research use only.

GENE DESCRIPTION	peptidyl-tRNA hydrolase 2
GENE ID	51651
ORIGIN SPECIES	Human
ACCESSIONS	BC006807
VECTOR	pDONR223
NCBI SYMBOL	PTRH2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccctccaaatccttggttatggaatattggtcatccagtcactcggcttggtgttggaattgcttggtcatgtgcctggcctggagccttcgagtatccttgggatctccccaaaagcaagcagcagcaagacacacagatactgaaagtgaagcaagcatcttggagaca
gcggggagtagaagatgattctgtggttcgaaatgactaaagatgggaaaagggaagtggtgctgccagtgctctcatgctgctgtttcagcctacaagcagattcaagaagaatctgaaatgctcaacaatgggaatactgtggccagcccaaggtggtgcaaaagctcctgat
gaagaaccctgattgcattattggcccatgcaaaaatgctgggactgactgtaagttaattcaagatgctggacgtactcagattgcaccaggctcctcaaaactgctagggaattggccaggaccagcagacctaattgacaaagtcactggtcacctaaacttactag