

SFT2D3-SFT2 domain containing 3 Gene

Reference: PTXBC006808

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	vesicle transport protein SFT2C; SFT2 domain-containing protein 3; SFT2 domain containing 3		

For research use only.

GENE DESCRIPTION	SFT2 domain containing 3
GENE ID	84826
ORIGIN SPECIES	Human
ACCESSIONS	BC006808
VECTOR	pDONR223
NCBI SYMBOL	SFT2D3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcttaacgggtctctcggaaatcctgcaaatagaaagataattctagatccggaatcctgtatctggtggaaccatggatttctacaagctgaattatcctcattgtatag