

RPP30-ribonuclease P/MRP 30kDa subunit Gene

Reference: PTXBC006991

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TSG15; ribonuclease P protein subunit p30; RNase P subunit 2; RNaseP protein p30; ribonuclease P/MRP 30kDa subunit; ribonuclease P/MRP subunit p30		

For research use only.

GENE DESCRIPTION	ribonuclease P/MRP 30kDa subunit
GENE ID	10556
ORIGIN SPECIES	Human
ACCESSIONS	BC006991
VECTOR	pDONR223
NCBI SYMBOL	RPP30
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgggtgttgcagattggacgtcgagcgggttctgacctaaggtctgcgcggactgtggagacagccgtcaccttggtctatcagttgtgtctatcatatcgttgacttaaggaaaagaaacaggaaattgaaaaccagtagctgtttctgaactcttcacaacttgcctaa
ttgtacagggaatacaagaccaattaaaatttaactagattaaacattattgtctcggatccatctcactgcaatgttttgagagcaactcttcaaggcccggtctatgatgtgtgtgcagttttccaagacagaaaagcttttcatattgttcacacatttagatgtgatttagctgcata
actgtaacagagaaactaccattttacttcaaaagacctctctaatatgtggcgattgaccgagggcctgttgaacttgcctatagccctgctatcaaaagactccacaatgagaaggatatacaatttccagtgccctcaatttgatgcaaatctgcaaggaaaagaaatgaattatctagtgct
gcagaaaggcctttagaataagaggccatagacgtggcaaatctaggtgtgttgggtctctgaaagtgcgccaaggctgcgggtgtcccaactgccgagcagcgtctccatggagaaactagaaaactgcttttgaattatctctacagtgaaagaaacctcgccatca
gaaggagatgaagattgtctccagcttccaagaaagccaagtgtgagggtga