

CRYAB-crystallin, alpha B Gene

Reference: PTXBC007008

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CMD1II; CRYA2; CTPP2; CTRCT16; HEL-S-101; HSPB5; MFM2; alpha-crystallin B chain; epididymis secretory protein Li 101; heat shock protein beta-5; heat-shock 20 kD like-protein; renal carcinoma antigen NY-REN-27; rosenthal fiber component; crystallin alpha		

For research use only.

GENE DESCRIPTION	crystallin, alpha B
GENE ID	1410
ORIGIN SPECIES	Human
ACCESSIONS	BC007008
VECTOR	pDONR223
NCBI SYMBOL	CRYAB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacatgccatccaccacctggatccgcccctcttcttcttcaactccccagccgctctttgaccagttctcgagagcacctgtggagctgatcttcccgacgtctacttccctgagtccttctacctcggccacctcttctcgggcaccagctggttgacac
tggactctcagagatgcgctggaaaaaggacaggtctctctgcaactggatgtgaagcacttctcccagaggaactcaagttaagggttgggagatgattgagtgcatggaaaacatgaagagcgcagagatgaacatggttcatctccagggagttccacaggaataccgg
atcccagctgatgtagacctctaccattacttcatccctgcatctgatgggtctcactgtgaatggaccaaggaacaggtctctggccctgagcgcaccattccatcaccctgaagagaagcctgcttcaccgcagcccccaagaatag