



Reference: PTXBC007009

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	carboxypeptidase A2; carboxypeptidase A2 (pancreatic)		
For research use only.			

GENE DESCRIPTION	carboxypeptidase A2 (pancreatic)
GENE ID	1358
ORIGIN SPECIES	Human
ACCESSIONS	BC007009
VECTOR	pDONR223
NCBI SYMBOL	CPA2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]