



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

FABP3-fatty acid binding protein 3, muscle and heart (mammary-derived growth inhibitor) Gene

Reference: PTXBC007021

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FABP11; H-FABP; M-FABP; MDGI; O-FABP; fatty acid-binding protein, heart; fatty acid binding protein 11; fatty acid binding protein 3, muscle and heart; heart-type fatty acid-binding protein; mammary-derived growth inhibitor; muscle fatty acid-binding prot		

For research use only.

GENE DESCRIPTION	fatty acid binding protein 3, muscle and heart (mammary-derived growth inhibitor)
GENE ID	2170
ORIGIN SPECIES	Human
ACCESSIONS	BC007021
VECTOR	pDONR223
NCBI SYMBOL	FABP3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggagacgtttctcgtggccactggaagcctatgtagacagcaagaatttcgatgactacatgaagtctactcgtgtgggttttctaccaggcaggtggccacatgaccaagcctaccacaatctcgtaaaagaatggggacattctcacctaaaaacacacagcacttcaagaaca
cagagatcattcttaagtgtgggttggaattcgtagagacacagacagatgacaggaaggctcaagtcattgtgacatgtagtggagggaacattgtcacttcagaaatgggacgggcaagagaccacactgtgctgggagcttaattgatgtaaaactctctgacactaccacag
gcactcgcctctccactccactctgaataaacaggcactga