

SAA4-serum amyloid A4, constitutive Gene

Reference: PTXBC007026

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C-SAA; CSAA; serum amyloid A-4 protein; constitutively expressed serum amyloid A protein; serum amyloid A4, constitutive		

For research use only.

GENE DESCRIPTION	serum amyloid A4, constitutive
GENE ID	6291
ORIGIN SPECIES	Human
ACCESSIONS	BC007026
VECTOR	pDONR223
NCBI SYMBOL	SAA4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggctttcacaggcattgtttctgctccttggtcatgggagtcaccagtgaagctggcgctcttttcaaggagctctccaaggggtgggacatgggcagagcctattggacataatgatatccaatcaccaaaattcaacagatatctctatgctcgggaaactatgatgct
gccccaaagaggacctgggggtgctggtgctgtaactcatcagccgttcagggtctatcttcagggttaataatagactactattttggaacagcagcactgtattggaggactcgaagtccaacgagaagctgagggaatggggccgagtggtgcaagaccccgaccgcttcaga
cctgacggcctgcctaagaataactga