

ZMYM5-zinc finger, MYM-type 5 Gene

Reference: PTXBC007048

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC050; MYM; ZNF198L1; ZNF237; zinc finger MYM-type protein 5; zinc finger protein 198-like 1; zinc finger protein 237; zinc finger, MYM-type 5; zinc finger MYM-type containing 5		

For research use only.

GENE DESCRIPTION	zinc finger, MYM-type 5
GENE ID	9205
ORIGIN SPECIES	Human
ACCESSIONS	BC007048
VECTOR	pDONR223
NCBI SYMBOL	ZMYM5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaaaatgttcagtgaggattagagtgactgaacagactcctgctttattagggaatatggccatggcaactagtcctcatggacataggggattcatttggtcaccagcttgctcttagtcagtagatctaggaactcaccagtgaagatgatgatgatgatgatgtgtgtttatt
gaatctatacaacctccttcaatttctgctccagcaatagctgatcaagaacttcattttgcatcatcaaaaaatgaaaagcctcaaggaaattattctgtaattcctcctcttcaagagatttggcatctcagaagggaataatagtgagacaattgttattgatgaagaggacgtagaa
acaaatggaggagcagagaaaaagtctctctttttatcgaatgggactctcctggaaactaaaaacaaacacgatttggatttctccacttcagctcttcaagaagtaagaccaagactggagtaagaccttttaaccctggtagaatgaatgtggcaggagacttattcagaatggagaa
tttgcactcatcatagctcctgagatgcattcaaaagaaggctaattgtcattctccagtag