

UBE2V2-ubiquitin-conjugating enzyme E2 variant 2 Gene

Reference: PTXBC007051

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DDVIT1; DDVit-1; EDAF-1; EDPF-1; EDPF1; MMS2; UEV-2; UEV2; ubiquitin-conjugating enzyme E2 variant 2; 1 alpha,25-dihydroxyvitamin D3-inducible; DDVit 1; MMS2 homolog; enterocyte differentiation promoting factor; enterocyte differentiation-associated facto		

For research use only.

GENE DESCRIPTION	ubiquitin-conjugating enzyme E2 variant 2
GENE ID	7336
ORIGIN SPECIES	Human
ACCESSIONS	BC007051
VECTOR	pDONR223
NCBI SYMBOL	UBE2V2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggctccacagaggttaaagttcctcgttaatttctgttgaagaactgaagaagacaaaaggagtaggcgacggtacagtagctggggccttgaagatgatgaagatatgacactacaaggtggacaggcatgattattggccaccaaggacaaattatgaaaacag
aatatatagcctgaaagtagaatgtggacctaataaccagaagctcctccgtcagtttagatttgaacaaaaaataatgaacggaataaataattccagtgggatgggtggatgcccgagcataccagtgtagcaaaatggcaaaatcatatagcattaaagtgtacttcaagagctaag
acgtctaagtagtccaagaaaatgaagctccacagccaccagaaggacaaacatacaacaattaa