

TBC1D7-TBC1 domain family, member 7 Gene

Reference: PTXBC007054

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MGCPH; PIG51; TBC7; TBC1 domain family member 7; TS complex subunit 3; cell migration-inducing protein 23		

For research use only.

GENE DESCRIPTION	TBC1 domain family, member 7
GENE ID	51256
ORIGIN SPECIES	Human
ACCESSIONS	BC007054
VECTOR	pDONR223
NCBI SYMBOL	TBC1D7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgactgaggactctcagagaaacttctcgtcagtatattatgagaaagtgggttctcgtggagttgaagaaaagaatacattagaaattctcctaaagatgaccgtctgatactgagaaacttctacttttagtcagaggttccctctccgtccatgtaccgtgcattggatggaggtgctt
ctaggaaatcttgcctccacaccacagagtcctcatgccaaggtgatgatgatcgttaaggagcagtagtggatgtcctcatgcccgtgaaagtcgttctgttgaagtgcacacacaggtgaggtctatctccgcatgtatcagctggagctctgggaagttacctcgaagtcctctttcc
actggagccagatgatgaagtgttcttccatagctaaagccatggaggaatggtggaagatagtgctgactgttactggatcacccgacgcttctgtaaccaatataatcacaagtaaccgggattcttgcctcagtgccaaagcgtttgaacaatactgaatctggaagatggcaga
ctgctgactcatctgaggaatgttccgcgccgccaactcttatgatctctggttcaagaggtgcttgcgggatgttgcctgaatccagtttacagagggttgggataaagttgtgagtggaatcctgtaagatcctagttttgtagctgtcgaatttttataacctttaaataaaagtattg
gcactgaacagtgagagaagatacaaaagtttctgaaaattatccaggacagctcagacgcgagctgtgagcaaggccattgactgtggcacaacacgttgggaccccggtccattcaagctga