

CLPS-colipase, pancreatic Gene

Reference: PTXBC007061

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	colipase; colipase, pancreatic; pancreatic colipase preproprotein		
For research use only.			

GENE DESCRIPTION	colipase, pancreatic
GENE ID	1208
ORIGIN SPECIES	Human
ACCESSIONS	BC007061
VECTOR	pDONR223
NCBI SYMBOL	CLPS
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagaagatcctgatcctcctgtgtgcccctctgtggcctatgcagctcctggccccgggggacattatcaacctggagaacggtgagctctgcatgaatagtgccagtgtaagagcaattgctgccagcattcaagtgcgctggcctggccccgtgcacatccatggccag
cgagaacagcgagtgctctgtcaagacgtctatgggatttactacaagtgtccctgtgagcgtggcctgacctgtgagggagacaagaccatcgtgggtccatcaccaacacaaacttggcatctgccatgacgctggacgtcctcaagcagtga