

RFK-riboflavin kinase Gene

Reference: PTXBC007069

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RIFK; riboflavin kinase; 0610038L10Rik; ATP:riboflavin 5'-phosphotransferase; flavokinase		

For research use only.

GENE DESCRIPTION	riboflavin kinase
GENE ID	55312
ORIGIN SPECIES	Human
ACCESSIONS	BC007069
VECTOR	pDONR223
NCBI SYMBOL	RFK
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccccgagcggactgcattatgaggcacctgccttactctgccgggtcaagtgtgcgggctcggccggtcccaagcagctgggcatccccacagctaatcttctgagcaagtggtagataatctccagctgatatccactggtatttactatgggtggccagtgttgaa
gtggagatgtccataagatgggtgagcataggatggaacccattacaagaatacgaagaagtctatggaacacatatcatgcatacctcaagaggacttctatgggaaatcctcaatgtcgccattgttgctacctgagaccagaaaagaacttgaattcttagagtcacttattc
agcaattcaagggtatattgaagaagctaagaacgactagagttaccagaacatttgaaaatcaagaagacaatttctccaggtttctaaagcaaaataatgaatggccactga