

ZMYM6-zinc finger, MYM-type 6 Gene

Reference: PTXBC007070

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Buster2; MYM; ZBED7; ZNF198L4; ZNF258; zinc finger MYM-type protein 6; transposon-derived Buster2 transposase-like protein; zinc finger protein 258; zinc finger, BED-type containing 7; zinc finger, MYM-type 6; zinc finger MYM-type containing 6		

For research use only.

GENE DESCRIPTION	zinc finger, MYM-type 6
GENE ID	9204
ORIGIN SPECIES	Human
ACCESSIONS	BC007070
VECTOR	pDONR223
NCBI SYMBOL	ZMYM6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaagaaccttggatggtgaatgtggcaaacgagtggtaccacagcaggagcttctggacaaaattaaagaaccagacaaatgctcaagagatggtatgtgtccacagccaaaactcaagaaagtaattgaaattgggtgtgtcttcagttatgagagacctattgccca
gcagttgaaccaggcttcagcttcttttgcacatctggcccaagtggtgtgtcttcagttccagctgtgtgtatgaaggtttttgtctgtgtgttaaaaaaatgctttataaggccaaactgcatatcataagacaggatctactcagctctctgtccacacgatgcaccagacattct
tcacctgctgctgccacctctcccaagaaaacctgcacaaactgctcgaagtataaaattcttaacatccctttttactttaccttttttag