

HOXB13-homeobox B13 Gene

Reference: PTXBC007092

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PSGD; homeobox protein Hox-B13; homeobox B13		
For research use only.			

GENE DESCRIPTION	homeobox B13
GENE ID	10481
ORIGIN SPECIES	Human
ACCESSIONS	BC007092
VECTOR	pDONR223
NCBI SYMBOL	HOXB13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagcccgcaattatgccaccttgatggagccaaggatcgaaggcttctggggagcgggagggggggaatctggtcgccactccctctgaccagccaccagcgccctacgctgatgctgtgcaactatgcccccttggatctgccagctcggggagccgc
caaaagcaatgccaccatgccccctgggggccccaggggagctccccagctccctgcttacttggaggcgggtactactcctgccagtgctccggagctcgtgaaacctgtgcccagcgagccacctggccgctaccccgccgagactccacggcggggaag
agtaccccgccgccccactgagtttgccttctatccgggatccgggaacctaccagcctatggccagttacctggagctgtctgtggtgcagactctgggtgctcctggagaaccgacatgactccctgttgcctgtggacagttaccagctctgggctcgtggtggaaca
gccagatgttggcaggagagaacagaccaccaggtcccttttggaggcagcatgtgcagactccagcgggcagcaccctcctgacgcctgcgccttctgctcgccgcccgaagaacgcatccgtacagcaaggggagcttgcgggagctggagcgggagtgccgctaaca
agttcatccaaggacaagagcgcaagatctggcagccaccagcctctggagcgccagattaccatctggttcagaaccgccgggtcaagagaagaaggttctcgccaaggtgaagaacagcgtacccttaa