

RQCD1-RCD1 required for cell differentiation1 homolog (S. pombe) Gene

Reference: PTXBC007102

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RQCD1; CAF40; CT129; RCD-1; RCD1; cell differentiation protein RCD1 homolog; RCD1 required for cell differentiation1 homolog; cancer/testis antigen 129; protein involved in sexual development; CCR4-NOT transcription complex subunit 9		

For research use only.

GENE DESCRIPTION	RCD1 required for cell differentiation1 homolog (S. pombe)
GENE ID	9125
ORIGIN SPECIES	Human
ACCESSIONS	BC007102
VECTOR	pDONR223
NCBI SYMBOL	RQCD1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcacagcctggcagggctgcgctgtgcctactacactggcacaagtggatagagaaaagatcactcagtgatcaatgagctgtccagctcctgagactagggaatcttctgctggagctaaagaaagcagaaatctgttctgacctgcacctatgctggcattcatttgg
actattgcagcactttacaggaattgtaaatattatccatctatcaacccacccacttgacagcacaccagtctaacagagtttcaatgctctggcattactgcaatgtgtagcatcacatccagaaccaggtcagcgtttctgcagcacatcccactttttgtacccttttgcaca
ctgtcagcaaacacgtcccttgagtatctccggctcaccagcctggagtactggggccctgggtaaacagatgaacaagaatcaacttttataaacacagaattatccctttatgtttgcgaattatggaatctggaagtgaacttctaaacagttgccacattcctccaga
agatctgttagatgacactgtttggcttatatatgtcagacgtatgagcgtttctccatgttgccatgctgggtaagatgtctcagctatccaagagcctctgccctctgctgaagcatgtagtgagatgttacctcgaacttcagataaccaggtttcagattgacttctct
ggctatcttcaagaaaatga