

C11orf58-chromosome 11 open reading frame 58 Gene

Reference: PTXBC007103

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IMAGE145052; SMAP; small acidic protein; chromosome 11 open reading frame 58		
For research use only.			

GENE DESCRIPTION	chromosome 11 open reading frame 58
GENE ID	10944
ORIGIN SPECIES	Human
ACCESSIONS	BC007103
VECTOR	pDONR223
NCBI SYMBOL	C11orf58
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtgtgtccagagagtgctcaccgcatgggtgaagcgttcagctccccagacgacgatctggatctagcaattgggaggcagcagacttgggtatgaagagagaaaaaagttcttgagacttatgggtgcaggaaagaaacatactggctgtttataggagatc
acaaatcaacatctcactccgaaccgggaagaagacaagaaaataatgaagaactggagtctcaatcagcaagatggacagtaaatatcaggaagatcggcgacattgtgacttggcttcagttaggtagaagaccatgatggagaagtgatgtggctggagatgatgat
gatgacgatgatgattcacctgatctgaaagtcagatgattctgaaagcattcagagtcagagaaagaagaatctgctgaagaactccaagctgctgagcacctgatgaagtggagatcccaaaacaaaaaagatgcaaaaagcaattataaatgatgtttgttaaatcagtggt
tataa