

ZNRF1-zinc and ring finger 1 Gene

Reference: PTXBC007235

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E3 ubiquitin-protein ligase ZNRF1; NIN283; nerve injury gene 283; nerve injury-induced gene 283 protein; zinc and ring finger 1, E3 ubiquitin protein ligase; zinc and ring finger protein 1; zinc/RING finger protein 1; zinc and ring finger 1		

For research use only.

GENE DESCRIPTION	zinc and ring finger 1
GENE ID	84937
ORIGIN SPECIES	Human
ACCESSIONS	BC007235
VECTOR	pDONR223
NCBI SYMBOL	ZNRF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggggcaagcagagcacggcgcccgctcccggggcccctcccggggtctccaccgatgacagcgccgtgccgcgcccgggagggggcggccatttcgggcactaccggacggcgccggggccatggggctgcgcagccgctcggtcagctcgggtgcaggcatg
ggcatggaccocagacggcgccgggggggtgccccttggcctctacaccccgctcccggggacccggcgactccgagaggcgccggcgggaggggtctgcgtccgactccactatgccatggcaatggtaccaggagacggcgccggtcaccatagagacgggat
gctgtacctgggctcccgagcctcgtggcggtgctctacctctgcacatcgacccaggtggttcagctcgcatagtggttcaagtgcccatgtgctcaagctctgtgcttctgacgagatggaatgcactttataatgtgttgagcaaacctgcctctctacaacgatgatgtgt
gactaaagacggggtagtgtgtgatctgcctggaggagctgtgcagggggacacgatagccagctgccctgcctgtgcatctatcacaagctgcatagactcgtggttgaagtgaacagatctgtccggaacacctgcggactga