



Reference: PTXBC007263

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NIZP1; ZFP496; ZKSCAN17; ZSCAN49; zinc finger protein 496; NSD1 (nuclear receptor binding SET-domain containing 1)-interacting zinc finger protein 1; zinc finger protein with KRAB and SCAN domains 17		

For research use only.

GENE DESCRIPTION	zinc finger protein 496
GENE ID	84838
ORIGIN SPECIES	Human
ACCESSIONS	BC007263
VECTOR	pDONR223
NCBI SYMBOL	ZNF496
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]