



Reference: PTXBC007287

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CR53; ZKSCAN21; ZSCAN53; zinc finger protein 213; zinc finger protein with KRAB and SCAN domains 21		

For research use only.

GENE DESCRIPTION	zinc finger protein 213
GENE ID	7760
ORIGIN SPECIES	Human
ACCESSIONS	BC007287
VECTOR	pDONR223
NCBI SYMBOL	ZNF213
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1