

ZNF414-zinc finger protein 414 Gene

Reference: PTXBC007310

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZFP414; zinc finger protein 414		

For research use only.

GENE DESCRIPTION	zinc finger protein 414
GENE ID	84330
ORIGIN SPECIES	Human
ACCESSIONS	BC007310
VECTOR	pDONR223
NCBI SYMBOL	ZNF414
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggagaacccctcaggcccatccagacatgctgccactgcagagcccagctccagtgagaccgacaaggaggtgtgtcccggtgtgtccagctgcagccctctctctccatgtcggaggagccagccctgagcaggcagccacaccgagtggtggaacgtgg
aggggctggaggatgcggcagggtctctcccccagcccagacagctgccagctggcccggaaccagccctggcctgaccagcatagctccgggaccagcgaggacctgcggcctcccagacgagccacctccagggaagcaaatccctgtctccagccctggctgtg
cctcagtttccagcgtcgtgacgtggcagcatctgcgaacccactgcccggccacacagtccttggaggcaagctcttccgctgtcagccctgagctgcaccgagaccttcccagcatgcaggagctggtggtcagagcaaaactgactacaaaccaatcgctactcaa
gtgtgagaactgcctcgtcgttccgcagcaccgctcgtcttcaagcacctgcatgtttgctgcggagcagcagagcccagcccgccgacccccggccctggagcagagagccgcccgcggcggagcggcccccggaggttgaccccgcatcagcgcgggctgccc
ttccgctgctggaaccttcacgaccccgccctgccccacggaccgttctgcccacttgaaccttgcccttggcctaagcccccgccctgcgcccctctgtccgctgcacccggcgccggtctccagcgccgctctgaaaaagagccaaggtgctgga
gcagcccccgaagaccagggcggtccgacgcgcccctcaggtgctgacgtga