



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

KIRREL2-kin of IRRE like 2 (Drosophila) Gene

Reference: PTXBC007312

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FILTRIN; NEPH3; NLG1; kin of IRRE-like protein 2; kin of irregular chiasm-like protein 2; nephrin-like gene 1; nephrin-like protein 3; kin of IRRE like 2 (Drosophila)		

For research use only.

GENE DESCRIPTION	kin of IRRE like 2 (Drosophila)
GENE ID	84063
ORIGIN SPECIES	Human
ACCESSIONS	BC007312
VECTOR	pDONR223
NCBI SYMBOL	KIRREL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr