

YKT6-YKT6 v-SNARE homolog (S. cerevisiae) Gene

Reference: PTXBC007319

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	YKT6 v-SNARE homolog (S. cerevisiae); YKT6, S. cerevisiae, homolog of; YKT6 v-SNARE protein; SNARE protein Ykt6; synaptobrevin homolog YKT6; R-SNARE		

For research use only.

GENE DESCRIPTION	YKT6 v-SNARE homolog (S. cerevisiae)
GENE ID	10652
ORIGIN SPECIES	Human
ACCESSIONS	BC007319
VECTOR	pDONR223
NCBI SYMBOL	YKT6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagctgtacagctcagctcctctacaaggcaggaaggtggtgctcctcaagccgcatacgtgtgtctcctcagcttttccagagatccagcgttcaggaattcatgaccttcacgagtcgaactgattgtggagcgtcatcgaaggcactagagcttctgcaagaac
aagactatctgtgccagctctacgtccggaatgatagctcttgaggtgtggtcattgctgacaatgaataccatcccgggtggcctttaccttgctggagaaagtactagatgaattctccaagcaagtcgacaggatagactggccagtaggagatccctgctacaatccattaccagccct
ggatggtcacctcagtagatagacagaaccacgagaagctgatccatgactaaagtgcaggccgaactagatgagacaaaatcattctgcacaacacatggagctctgttagagcgaggtgagaagctagatgactgggttccaatccgaggtgctgggaacacagctaaagc
cttctataaaactgcccgaacaaaactcatgctgtgccatcatgtga