

ANXA10-annexin A10 Gene

Reference: PTXBC007320

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ANX14; annexin 14; annexin-10		

For research use only.

GENE DESCRIPTION	annexin A10
GENE ID	11199
ORIGIN SPECIES	Human
ACCESSIONS	BC007320
VECTOR	pDONR223
NCBI SYMBOL	ANXA10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttttggagactatgtgcaaggaaacatctccagctcccaattcaatccataatggatgcccaatgctaggaggagcactcaaggattgactgtgacaaagacatgctgatcaacattctgactcagcgtgcaatgcacaaaggatgattgcagaggcataccagagca
tgtatggccgggacatgattgggatatgaggagcagctttcggatcactcaaaagatgtgatggctggcctcatgtaccaccaccactgtatgatgctcatgagctctggcatgccatgaaggaggtaggcactgatgagaattgcctcattgaaatactagctcaagaacaaatggag
aaattttcagatgcgagaagcctactgcttcaatacagcaataacctcaagaggacattttcagagacctcaggacacttcagagatactctcatgaacttggtccaggggaccagagaggaagatatacagacctgcgatggctgctcaggatgcaatggtcctatgggaagcc
tgtcagcagaagacggggagcacaaaacctgctgcaaatgatcctgtgcaacaagagctaccagcagctcggctggtttccaggaaattcaaaatattctgggcaagatatggttagatccattaatgaatgttatgatggatacttcaggagctgctggttgcaattgtctctgtgttc
gagacaaccagcctatgttcttatgattatagtgcaattcatgcttgggttcataataaaactgtaacagcaggttcattgccagaagtgaatagacctgctgaccataaggaaacgatacaaaagagcgatatggaaaatccctattcatgatacagaaatttgcctcaggccatta
taagaagcactgcttccatctgtgctggtgatgctgaggactactaa