

NDUFC2-NADH dehydrogenase (ubiquinone) 1, subcomplex unknown, 2, 14.5kDa Gene

Reference: PTXBC007323

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	B14.5b; CI-B14.5b; HLC-1; NADHDH2; NADH dehydrogenase [ubiquinone] 1 subunit C2; NADH dehydrogenase (ubiquinone) 1, subcomplex unknown, 2, 14.5kDa; NADH-ubiquinone oxidoreductase subunit B14.5b; complex I subunit B14.5b; complex I-B14.5b; human lung cancer		

For research use only.

GENE DESCRIPTION	NADH dehydrogenase (ubiquinone) 1, subcomplex unknown, 2, 14.5kDa
GENE ID	4718
ORIGIN SPECIES	Human
ACCESSIONS	BC007323
VECTOR	pDONR223
NCBI SYMBOL	NDUFC2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgcacggcggaaccagaaacccctacggtttctgcccgatgagggccggagcctgccccgcccaagctgaccgacccggcctcctacatcggtcttgggctactgctccggcctgattgataacctgatccggcgaggccgatcgacggctggttgcacgcca
gcttctatataftacggccttttttctgctggaattatctgtgaaacgtgaagactacctgtatgctgtgaggacccgtgaaatgtttggatatafgaaattaccagaggatttctgaagaagataagaaaacatatggtgaaattttgaaaaattccatccaatcgttga