

BYSL-bystin-like Gene

Reference: PTXBC007340

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	by the ribosomal protein s6 gene, drosophila, homolog-like; bystin like		

For research use only.

GENE DESCRIPTION	bystin-like
GENE ID	705
ORIGIN SPECIES	Human
ACCESSIONS	BC007340
VECTOR	pDONR223
NCBI SYMBOL	BYSL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgacagcagcggccatcatcagagggtgtgtggaccctgaggatgagcgtgccatagatgttcatgaacaagaaccctctgccaggcaccctgctgacatcatatggagaagctgactgagaagcagacagaggtgagacagtcagaggtcggcttcc
ctatgccccagctggacccccgggtcctagaaggtgtacaggggtccgggaggtattatctaagaccgcagtggaaaactgcccaaggcatttaagatcatcctgcacttccaactgggatcaaatcctctacgacagagccggagcctggactgcagctccatgtaccag
ccaccaggattttgctcctaactgaaggaaacgatggcccagcgtcttacaacctgtctctccctcagtagagatgacgttctgaatacaaacgactcaactcactctacatggctctcaagaaggccctttcaaacctggagcctggttcaaaaggatcctgattccactgt
gcgagctggcactgtaccctccgggaagccatcattgtgggtagcatalcaccaggtgtccatccctgtgtgctcactccagtcggccatgctgaaaatgctgagatggaalacagcgggtccaacagcatctccgactgctgtgataagaagtatgcactgcctaccgggt
gctggatgccctagcttccactcctgggttccggacagagaagcgtgaactcctgtgtgtggcaccagtgccctctgacttgggtccagcgtacaaggccgacttggccacagaccagaagaggccctctagaactgctccggctgcagccccatccacagctatgccgaa
atcaggcgtgagctcagagtcagctcccccgatgtggaagatgttccatcaccgtggagtga
```