

DYDC2-DPY30 domain containing 2 Gene

Reference: PTXBC007374

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DPY30 domain-containing protein 2; DPY30 domain containing 2		

For research use only.

GENE DESCRIPTION	DPY30 domain containing 2
GENE ID	84332
ORIGIN SPECIES	Human
ACCESSIONS	BC007374
VECTOR	pDONR223
NCBI SYMBOL	DYDC2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaactaactacgtgaagaggtgcttggaaattgcctggccagggcactggcagaggtggcgaaggttcggccagtgacccaatagaatacctggctcactggctttatcattacaggaaaacagcaaaagcaaaagaagaatagggaagaagatccactgcaggagg
aatatgacagtagcctcaaggaaatggaaatgacagaatgctgaacagggaagatcagattcaacagaactgtgaaaagtgtcacaagggaactgactctgaaactgtttccaggaagaagaccatattcatgcaggaggacacaaaccccttgagaaggaggccttgaagcagg
aattcctgccaggtacttcagctgattcaggaatgcctcaacaggttcctctcagagctctgtgccagattgaccagaacttcaaaatgcccaagaataaattacaaggaggtttcagcatgaagttgctcatgaaatgcctcctggtccaaatctcttttag