

FBXW4-F-box and WD repeat domain containing 4 Gene

Reference: PTXBC007380

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DAC; FBW4; FBWD4; SHFM3; SHSF3; F-box/WD repeat-containing protein 4; F-box and WD-40 domain protein 4; F-box and WD-40 domain-containing protein 4; F-box/WD repeat protein 4; dactylin; F-box and WD repeat domain containing 4		

For research use only.

GENE DESCRIPTION	F-box and WD repeat domain containing 4
GENE ID	6468
ORIGIN SPECIES	Human
ACCESSIONS	BC007380
VECTOR	pDONR223
NCBI SYMBOL	FBXW4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccctggatgcagctagaggatgattctctgtacatacccaggctaatttcacctcgccctaccagttccgtccagatggtgccagctgaatcgtcgccctctgggagcttttctgggcatgaggacgtttgccactttgtctggccaactcgcattattgttagtcaggaggggat
gggaagattggcattcataagattcacagcacttcactgcaagtactcggctcatgaacaggaggtgaactgtgtggattgcaagggggcatcattgtgagtggtccaggacaggacggccaaggtgtggcctttggcctcaggccggctgggagctgttacacaccatccag
actgaagaccgagctgtgtccattgctatcagccattactcagctctttgtgacaggggcgtgtgtggcgcactctcaccctgagaatctgggacctcaacagtgggcagctgatgacacacttgggcagtgacttccccaggggctggggtgtgtgatgtatgatgcc
ctttcacactgtctgtggctatgacacctatgttcgctactgggacctccgaccagcgtccggaatgtgtcagtgagtgaggaggagcccccagacagcaccctgtactgcctgcagacagatggcaaccacgtctggccacaggttctcctactacggtgtgtacggcgtgg
gaccggcgtcaaggcctgctgcacgcctccgctgacgtcgactccccctcagcagccgtgtgtactgcctgcgtctcaccaccaagcatctctatgctgccctgtcttacaacctccagtcctggattttcaaaacccatga