

ZCCHC17-zinc finger, CCHC domain containing 17 Gene

Reference: PTXBC007446

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC251; PS1D; nucleolar protein of 40 kDa; nucleolar protein 40; pnn-interacting nucleolar protein; zinc finger CCHC domain-containing protein 17; zinc finger, CCHC domain containing 17; zinc finger CCHC-type containing 17		

For research use only.

GENE DESCRIPTION	zinc finger, CCHC domain containing 17
GENE ID	51538
ORIGIN SPECIES	Human
ACCESSIONS	BC007446
VECTOR	pDONR223
NCBI SYMBOL	ZCCHC17
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaattcaggaaggcctgagaccatgaaaactgcctgctctacactattttcaaggagaggttgctatggtgacagactatgggaccttatcaaatcccaggctgtcgaagcaaggctgtggtccatcgaactcatatgtcatcctgtcgggtgataagccctctgagatagtaga
tgttgagataaagtgtggtgaagcttattgcccagagatgaaaaatgataataaaatgacccctccatgaaggttgcaatcaaggagctgggaaagaccttgatcccaacaatgttatcattgagcaagaagagagggcggagggcgtccaggattacactgggcagaaga
tcaccttgaggctgtctgaacactacctgcaagaagtgtggtgtaaaaggccactttgcaaaagattgtttcatgcaaccagggtgggactaaatactctctgatacctgatgaggaaaggaaaaggaaaggcaaaagtcagcagagtttgagaagcctgacctacaaggaatcctcta
gaaaaagaaagaggaagaagaaaaaagaacatagagataggaagtcatctgactctgacagctcagactctgagagtgatacaggcaaggcgaagcacatcaaaagacagcaggcagcaaaagaagaagaaaaagaagaagacacagaagaagcacaagga
gtga